








OFTHALMISCHE CHIRURGIE

Niet - resorbearbare hechtingsdraad

POLYPROPYLENE *BLUE*

monofilament

uncoated



Naald	Draaddikte	Lengte	Kleur	Order nummer	Stuks/ds
1/2 Circle - side cutting spatulla 2 x 6,2 mm - 200 µm	 10/0	30 cm	blauw	OPP106	12
3/8 Circle - side cutting spatulla 2 x 6,4 mm - 200 µm 2 x 6,4 mm - 200 µm 6 mm - 200 µm	 10/0 9/0 7/0	30 cm 30 cm 45 cm	blauw blauw blauw	OPP105 OPP95 OPP71 - 1	12 12 12
1/4 Circle 16 mm - 150 µm 2 x 8 mm - 300 µm	 10/0 6/0	10 cm 45 cm	blauw blauw	OPP103 OPP62	12 12
1/4 Circle loop 2 x loop 16 mm - 150 µm 2 x loop 15,3 mm - 150 µm	 10/0 10/0	20 cm loop 20 cm loop	blauw blauw	OPP103L - 2 OPP113L - 2	12 12
Straight 2 x 16,5 mm - 150 µm 2 x 16,5 mm - 150 µm	 10/0 9/0	20 cm 20 cm	blauw blauw	OPP101 OPP91	12 12
Straight loop 16,5 mm - 150 µm	 10/0	15 cm loop	blauw	OPP110	12
Straight needle + 3/8 circle needle 1 x 16 mm straight needle - 150 µm 1 x 4 mm needle 3/8 - 100 µm	 10/0	10 cm	blauw	OPP102	12



Niet - resorbearbare hechtingsdraad

DACLON *NYLON BLACK*

monofilament

uncoated

Naald	Draaddikte	Lengte	Kleur	Order nummer	Stuks/ds
1/4 Circle - Round bodiet - taper point 2 x 8 mm - 220 µm	 8/0	45 cm	zwart	ON108	12
1/2 Circle - side cutting spatulla 2 x 5,5 mm - 150 µm 1 x 5,5 mm - 150 µm single armed 2 x 6,2 mm - 150 µm 2 x 6,4 mm - 220 µm	 10/0 10/0 10/0 10/0	30 cm 30 cm 30 cm 30 cm	zwart zwart zwart zwart	ON108 ON1081 ON102 ON106	12 12 12 12

2 x 8 mm - 300 µm	6/0	45 cm	zwart	ON62	12
3/8 Circle 					
2 x 6,2 mm - 150 µm	11/0	20 cm	zwart	ON111	12
2 x 6 mm - 220 µm	10/0	30 cm	zwart	ON103	12
2 x 6,2 mm - 150 µm	10/0	30 cm	zwart	ON101	12
1 x 6,2 mm - 150 µm single armed	10/0	15 cm	zwart	ON1011	12
2 x 6,4 mm - 220 µm	10/0	30 cm	zwart	ON105	12
2 x 6,4 mm - 220 µm	10/0	30 cm	zwart	ON105 - 38	12
2 x 6,4 mm - 220 µm	10/0	45 cm	zwart	ON105 - 45	12
2 x 6,2 mm - 150 µm	9/0	30 cm	zwart	ON95	12
2 x 6,2 mm - 150 µm	9/0	45 cm	zwart	ON95 - 45	12
2 x 6,2 mm - 150 µm	8/0	30 cm	zwart	ON85	12
2 x 6,2 mm - 150 µm	8/0	45 cm	zwart	ON85 - 45	12
2 x 6,2 mm - 200 µm	7/0	30 cm	zwart	ON75	12
2 x 8 mm - 300 µm	6/0	45 cm	zwart	ON61	12
1/4 Circle 					
2 x 8 mm - 350 µm	5/0	45 cm	zwart	ON51	12

Niet - resorbeerbare hechtingsdraad

POLYESTER BRAIDED WHITE


multifilament *gedraaid* coating *silicone*

Naald	Draaddikte	Lengte	Kleur	Order nummer	Stuks/ds
1/2 Circle - side cutting spatulla 					
2 x 8 mm - 350 µm	5/0	45 cm	wit	OP52	12
3/8 Circle 					
2 x 8 mm - 350 µm	5/0	45 cm	wit	OP53	12
1/4 Circle 					
2 x 8 mm - 250 µm	7/0	45 cm	wit	OP71	12
2 x 8 mm - 300 µm	6/0	45 cm	wit	OP61	12
2 x 8 mm - 350 µm	5/0	45 cm	wit	OP51	12
2 x 8 mm - 350 µm	4/0	45 cm	wit	OP41	12

Niet - resorbeerbare hechtingsdraad

SILK VIRGIN BLUE

multifilament *gedraaid* uncoated

Naald	Draaddikte	Lengte	Kleur	Order nummer	Stuks/ds
3/8 Circle - side cutting spatulla 					
2 x 6,5 mm - 150 µm	8/0	45 cm	blauw	OS82	12
2 x 6,5 mm - 200 µm	8/0	45 cm	blauw	OS81	12

SILK BRAIDED BLACK

multifilament *gevlachten* coating *silicone*

3/8 Circle - side cutting spatula



2 x 6,5 mm - 200 µm	8/0	45 cm	zwart	OS81B	12
2 x 6,5 mm - 240 µm	7/0	45 cm	zwart	OS71	12
2 x 6,5 mm - 240 µm	6/0	45 cm	zwart	OS61	12
2 x 8 mm - 350 µm	5/0	30 cm	zwart	OS52	12

3/8 Circle - reverse cutting



2 x 12 mm	7/0	45 cm	zwart	8051512 - 2	12
2 x 12 mm	6/0	45 cm	zwart	8071512 - 2	12

1/4 Circle

2 x 8 mm - 250 µm	6/0	45 cm	zwart	OS62	12
2 x 10 mm - 350 µm	5/0	45 cm	zwart	OS51	12

